

**FOREVER** YOUNG


 **DR. MÜLLER**  
FOREVER YOUNG



Aviv International Pte Ltd | 06-10, CT Hub, 2 Kallang Avenue,  
Singapore 339407. | Tel: 67872848 (AVIV) | [www.avivwellness.com](http://www.avivwellness.com)

*Dr. Müller is a registered trademark of Ultrasun International*

DR.MÜLLER INFINITY  
PRODUCT LEAFLET 2018



# DR.MÜLLER INFINITY

## DIMENSIONS

Closed



Open



When placing the Dr. Müller Infinity make sure that the device is surrounded by 10 -15 centimetres of free space for air circulation.

## LAMP SPECIFICATIONS

Top lamps	16 x Dr. Müller Collagenic Light 120W E-Tronic   190 cm
Base lamps	16 x Dr. Müller Collagenic Light 120W E-Tronic   190 cm
Recommended exposure time	20 minutes


## TECHNICAL SPECIFICATIONS

Power supply	3 x 400V + PE   50/60 Hz
Plug and cable	16A CEE 5 x 2,5 mm <sup>2</sup>
Power	4.200W
Breaker/fuse (slow)	3 x 16A
Minimum net weight <sup>(1)</sup>	180 kg
Minimum crated weight <sup>(1)</sup>	215 kg
Boxed dimensions	227 x 114 x 148 cm
Special boxed dimensions	227 x 114 x 133 cm


## STANDARD EQUIPPED WITH



## ADDITIONAL OPTIONS



## CONTROL SYSTEMS



## AVAILABLE COLOUR

Dream White Metallic

<sup>(1)</sup> Minimum weight is based on a Dr. Müller Infinity with standard configuration.

Machine specific weight information available upon request.

<sup>(2)</sup> Options only available in combination with Music Pro.